SALVAFIX Slim

ENGLISH




SALVAFIX Slim Abdominal belt *for seated patients. Indicated to prevent injuries and falls from the chair. It consists of a belt with adjustable fitting by eyelets and a base for fixing the set to the bed. Both parts are joined by reinforced seams without articulation. Available with sewn perineal to prevent sliding down of the patient as well an avoid the risk of suffocation.*

Made with polyester webbing, inner reinforcement and stainless steel eyelets for extra strength. Soft-touch polyester trim for a comfortable, chafe-free fit. Safety latches made up of a stainless steel pin and a mechanism with ABS/PC casing. Magnetic key with PA housing and neodymium magnet. The mechanism is locked by inserting the pin through its hole and unlocked by placing the magnetic key on its upper face.

SEAT WAIST BELT

FEATURES



Soft to the touch edges



Colour code for sizes
M (Blue) - L (Green)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Slim

PRODUCT NAME: Seat waist belt

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the chair

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Code	Measures (cm)		
	Size	Contour Min./Max.	Width
SVF3060P	M	68 / 108	8
SVF3062P	L	74 / 114	8

COLOUR

- Natural white
- Grey

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 1 - Salvafix Slim - Seat waist belt
- 2 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE