SALVAFIX Pro



SALVAFIX Memory PRO perineal waist belt *with safety locks and magnetic key. Prevents patient self-injuries and falling. Includes sewn perineal band to prevent patient from sliding up and down.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch blue edging for a comfortable fit.
Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

PERINEAL WAIST BELT

FEATURES



Patented Memory®
Alphanumeric system



Adaptable to 80-100cm beds



Sewn sidebands



Sewn perineal band



Soft to the touch edges



Reinforced articulation



Double eyelets for XL sizes



Colour code for sizes
S (Yellow) - M (Blue)
L (Green) - XL (Red)



Bariatric model
XXL (Brown)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Memory Pro

PRODUCT NAME: Perineal waist belts

PATENT: EP2384726B1

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Measures (cm)			
		Waist Belt		Perineal band	
		Contour Min. / Max.	Width	Length	Width
SVF3001P	S	46 / 76	16	60	5
SVF3000P	M	68 / 108	18	80	5
SVF3002P	L	95 / 135	18	100	5
SVF3003P	XL	120 / 150	18	120	5
SVF3004P	XXL	135 / 200	27	120 (x2)	5

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING

- 1 - Salvafix Memory PRO - Perineal waist belt
- 3 - Safety locks (4 safety locks for L - XL, 8 u for XXL)
- 1 - Magnetic key
- 1 - Safety and instructions manual





MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

UK CA CE