SALVAFIX Standard

ENGLISH



SALVAFIX Standard Complete System. *The complete Salvafix system is the comprehensive solution for bedridden patient containment. It provides the highest level of safety for both patients and medical staff. The system consists of all necessary elements from the Standard range for total patient containment. These devices can also be used individually if needed.*

COMPLETE SYSTEM

COMPONENTS


Harness


Wristbands


Sewn perineal band


Waist Belt


Ankle bands


Safety locks
Magnetic keys

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Standard

PRODUCT NAME: Complete system

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

Basic UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size
SVF3121P	S
SVF3120P	M
SVF3122P	L

COLOUR

- Natural white
- Light blue
- Brown

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low
- Do not iron

PACKAGING

Grouped into a cardboard box:

- 1 harness
- 1 bed strap for harness
- 1 waist belt
- 2 wristbands
- 1 set of ankle bands
- 12 locks (13 for L)
- 4 magnetic keys
- Instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice