SALVAFIX Pro

ENGLISH



***SALVAFIX Memory PRO Complete system**. The complete set is the global solution of secure restraint for bedridden patients. It permits having the maximum level of security of the patient and the medical staff, because of bed restraint is a not easy issue. As consists of individual elements can be used as required at the time.*

All grouped in a sturdy carrying bag, capable of holding the entire system. It features a thermoplastic handle for a better grip and a closure with white injected zipper. It also includes tags for each size to be placed on the side of the bag, indicating the size of its contents.

COMPLETE SYSTEM

COMPONENTS


Harness


Wristbands


Sewn perineal band


Transportation bag


Waist Belt


Ankle bands


Safety locks
Magnetic keys



Colour code for sizes
S (Yellow) - M (Blue)
L (Green) - XL (Red)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Memory Pro

PRODUCT NAME: Complete system

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

Basic UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size
SVF3151P	S
SVF3150P	M
SVF3152P	L
SVF3153P	XL

COLOUR

- Natural white
- Light blue
- Brown

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low
- Do not iron

PACKAGING

Grouped in a carrying bag and fitted into a cardboard box:

- 1 harness with thorax band
- 1 bed strap for harness
- 1 waist belt
- 2 wristbands
- 2 ankle bands
- 1 bed strap for ankle bands
- 15 locks (16 for L and XL)
- 4 magnetic keys
- 1 carry bag
- 1 extendable key-ring
- Application stickers
- Instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE