SALVAFIX Standard



SALVAFIX Standard wristbands *allow to ensure the patient's hands to the bed. These are fixed to the bed structure by a long adjustable band.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch piping for a comfortable fit. Welded stainless steel rings. Comfortable padded wristbands with 5mm foam and 100% polyester soft fabric and hook/loop fasteners allows faster and tight application, increasing system security. Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

WRISTBANDS

FEATURES



Comfortable padded wristbands



Soft to the touch edges



Colour code for sizes
S (Yellow) - M (Blue) - L (Green)



Colour identifies limb position:
Brown = GROUND (Lower parts)
Light blue = SKY (Upper parts)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Standard

PRODUCT NAME: Wristbands

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Measures (cm)	
		Contour Min. / Max.	Width
SVF3221	S	9 / 15	8
SVF3220	M	14 / 20	8
SVF3222	L	19 / 25	8

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 2 - Salvafix Standard - Wristbands
- 2 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice