
# SALVAFIX Standard



***SALVAFIX Standard Wrist Bands Monoblock*** *allow to ensure the patient's feet to the bed. The set composes a strap fixed to the bed structure and two sewn wrist bands.*

*Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch piping for a comfortable fit. Welded stainless steel rings. Comfortable padded wrist bands with 5mm foam and 100% polyester soft fabric and hook/loop fasteners allows faster and tight application, increasing system security. Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.*

## WRIST BANDS MONOBLOCK

## FEATURES



**Comfortable padded wrist bands**



**Soft to the touch edges**




**Colour code for sizes**
S (Yellow) - M (Blue) - L (Green)



**Colour identifies limb position:**
**Brown = GROUND (Lower parts)**
**Light blue = SKY (Upper parts)**

## SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Standard

PRODUCT NAME: Wrist Bands Monoblock

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

BASIC UDI-DI: 843603535SVFNC

### SIZE CHART

| Ref. | Size | Measures (cm) | |
|---|---|---|---|
| | | Contour Min. / Max. | Width |
| SVF3231 | S | 9 / 15 | 8 |
| SVF3230 | M | 14 / 20 | 8 |
| SVF3232 | L | 19 / 25 | 8 |

### COLOUR

- Natural white
- Light blue

### COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
*LATEX FREE*

### LAUNDRY INSTRUCTIONS




- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

### PACKAGING




- 1 - Salvafix Standard - Wrist Bands Monoblock
- 4 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice