SALVAFIX Standard

ENGLISH




SALVAFIX Standard mitt bed *allow to ensure the patient's hands to the bed and prevents self-injuries and equipment manipulation. It is fixed to the bed structure by a long adjustable band and safety locks.*

Straps made of 100% polyester with inner reinforcement for a higher tensile strength. Stainless steel rings and eyelets. Mitt is made of comfortable 8mm foam with polyamide velour, viscose fabric and polyester breathable mesh fabric. Polyester soft to the touch edge all around the mitten for a comfortable fit. Hook/loop fasteners for a faster application. Fiber filled for limiting wrist movements.
Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

BED MITT

FEATURES



Comfortable padded interior and breathable mesh fabric



Mitt available in 2 sizes: M - L



Bed restraint band with safety locks



Soft to the touch edges

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Standard

PRODUCT NAME: Bed Mitt

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Units			
		Mitt	Band	Locks	Key
SVF3210	M	1	1	2	1
SVF3212	L	1	1	2	1
SVF3290	M	10	-	-	-
SVF3292	L	10	-	-	-
SVF3295	🛏	-	10	-	-

Mitt	
Size	Measures
M	14 / 25 cm
L	14 / 25 cm

Restraint band	
Size	Contour Min./Max.
🛏	14 / 25 cm

COLOUR

○ Natural white ● Light blue

COMPOSITION

Restraint band:
Polyester webbing
Stainless steel eyelets
Stainless steel rings

Mitten:
65% PES - 30% PA - 5% Viscose
Silicone hollow polyester fiber

WARNING

- Read and understand the safety and instructions manual before using Salvafix products. Failure to follow instructions and safety recommendations could cause serious injuries.

LAUNDRY INSTRUCTIONS

Restraint band:



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

Mitten:



- Do not wash
- Product for a single patient

PACKAGING

PACK
- 1 - Bed Mitt
- 1 - Restraint band
- 2 - Safety lock
- 1 - Magnetic key
- 1 - Instructions manual

SVF3210-12 — BOX 5 PACKS

SVF3290-92 — BOX 5 MITTS

SVF3295 — BOX 20 BANDS

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE