SALVAFIX Pro



SALVAFIX Memory PRO independent thorax harness *prevents the patient from slipping vertically up. A band around the thorax prevents from slipping and loosen the harness. Stops a violent or immediate incorporation avoiding the patient to injure themselves or others. This harness can be placed without being connected to a waist belt and also offers the possibility to place the harness straps in parallel position (uncrossed) to avoid touching the neck or avoid pressure on thesternum.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch piping for a comfortable fit. Safety locks are made of a Stainless Steel pin and a PC/ABS casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

INDEPENDENT THORAX HARNESS

FEATURES



Adaptable to 80-100cm beds



Soft to the touch edges



8cm width thorax band



Colour code for sizes
S (Yellow) - M/L(Blue) - XL (Red)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Memory Pro

PRODUCT NAME: Independent thorax harness

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

Basic UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Measures (cm)	
		Contour Min. / Max.	Width
SVF3441	S	40 / 80	8
SVF3440	M/L	75 / 115	8
SVF3443	XL	105 / 145	8

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING

- 1 - Salvafix Memory PRO - Independent thorax harness
- 6 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual



MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice