SALVAFIX Standard

ENGLISH




SALVAFIX Standard universal DUO *is a versatile and universal size restraint valid for use in wrists and/or ankles.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch piping for a comfortable fit. Welded stainless steel rings. Comfortable padded wristbands with 5mm foam and 100% polyester soft fabric and hook/loop fasteners allows faster and tight application, increasing system security. Safety locks are made of a Stainless Steel pin and a PC/ABS casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

UNIVERSAL DUO

FEATURES



Comfortable padded insoles with a soft touch and quick fastening using a doble Velcro® system



Soft to the touch edges



One size
Allows coverage from T/S to T/L



For upper and lower limbs

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Standard

PRODUCT NAME: Universal DUO

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Measures (cm)	
		Contour Min. / Max.	Width
SVF3800	One size	10 / 34	8

COLOUR

- Natural white
- Grey

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 2 - Salvafix Standard - Universal DUO
- 4 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE