MEDICOVERS absor

ENGLISH



Medicovers Absor ECO *ultra-absorbent underpads protect the mattress and bedding from liquids, while being breathable, thus keeping the patient comfortably dry. They are washable -resists 150 washes at 90 degrees- and reusable, assuming a great advantage over the disposable underpads because are more respectful with the environment and becoming more economical.*

Advanced Absor technology: Made of 4 layers of technical fabrics. First layer quickly transfers liquids to the lower layers keeping the top dry in a few minutes. The next layers retain liquids keeping the mattress protected and dry.

ABSOR ECO UNDER PADS

FEATURES

- Advanced absor technology
- Impermeable layers
- Reusable (washable up to 150 times)
- Breathable
- Hygienic
- Inexpensive

SPECIFICATIONS

PRODUCT FAMILY: MEDICOVERS Absor

PRODUCT NAME: Absor Eco Underpad

INTENDED USE: Protection for bed and mattress against liquids and body fluids.

BASIC UDI-DI: 843603535ABSOR9Y

SIZE CHART

Code	Measures (cm)		
	Width	Length	Wings
ABSORECO7090A	70	90	40

COLOUR

Green / Light blue White

COMPOSITION

Under pad: 83% Polyester - 10% Viscose - 7% Polyurethane
Wings: 50% Cotton - 50% Polyester

LAUNDRY INSTRUCTIONS



- Do not use fabric softener, peroxides or compounds with chlorine

- Machine wash 90ºC
- Do not bleach
- Tumble dry low
- Do not iron

Do not dry clean

PACKAGING

- **1 Pack:** 10 Medicovers Absor ECO underpads
- **1 Box:** 4 Packs (10 Medicovers Absor ECO underpads)



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE