MEDICOVERS absor

ENGLISH



Medicovers Absor PLUS *ultra-absorbent underpads are designed to protect the mattress and the bed's clothes, maintaining the patient mobility and comfortably dry at the same time. Made with 5 layers of technical fabrics sewn with wavy stitching to provide more durability.*


Includes a 40 cm side wings that can be easily adjusted underneath the mattress to obtain a better and prolonged position of the under pad.

Advanced Absor technology: Made of 5 layers of technical fabrics. First layer quickly transfers liquids to the lower layers keeping the top dry in a few minutes. The next layers retain liquids keeping the mattress protected and dry.

ABSOR PLUS UNDER PADS

FEATURES

- Advanced absor technology (5 Layers)



1- ALWAYS DRY - 100% PES
2 - HIGH DISPERSION - 100% PES
3- HIGH ABSORPTION - 80% PES / 20% Viscose
4 - BREATHABLE/IMPERMEABLE - 100% PU
5 - PROTECTION - 100% PES

- Absorption speed: First drop is immediately absorbed at less than 60mbar of pressure for a new product and improved up to a 27mbar pressure after 200 washing cycles according standard (UNE EN 20811).
- Absorption capability: Retains more than 3,10 L/m2 even after 200 washing cycles.
- Reusable (washable up to 200 times)
- Breathable
- Hygienic
- Inexpensive

SPECIFICATIONS

PRODUCT FAMILY: MEDICOVERS Absor

PRODUCT NAME: Absor PLUS Under pad

INTENDED USE: Protection for bed and mattress against liquids and body fluids.

BASIC UDI-DI: 843603535ABSOR9Y

SIZE CHART

Code	Measures (cm)		
	Width	Length	Wings
ABSORPLUS7090A	70	90	40

COLOUR

Green / Light blue White

COMPOSITION

Under pad: 86% Polyester - 8% Viscose - 6% Polyurethane
Wings: 50% Cotton - 50% Polyester

LAUNDRY INSTRUCTIONS




- Do not use fabric softener, peroxides or compounds with chlorine

- Machine wash 90ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING

- **1 Pack:** 10 Medicovers - Absor PLUS underpads
- **1 Box:** 4 Packs (10 Medicovers - Absor PLUS underpads)



MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE