BODYFIX

ENGLISH



Bodyfix pyjama cover *designed to provide a comfortable and very light restraint in passive, disoriented or with neurological affections patients, preventing them from accessing the diaper, tubes, probes, etc.*

It's used similarly to a sheet, with a versatile elastic adjustment under the mattress. It allows postural changes and let the patient sit and move on bed freely and comfortably. Longitudinal zipper with easy access to the patient for hygiene and release.

Two small side holes with zips allow the entrance of tubes, probes, catheters, cables, etc. Made in two types of fabrics, blue terry cloth for cold weathers and white piqué for warm environments.

PYJAMA COVER

FEATURES



2 adjustable elastic bands



Zip from neck to feet



2 lateral holes for medical devices



2 fabric options (WARM / COLD)

SPECIFICATIONS

PRODUCT FAMILY: Bodyfix

PRODUCT NAME: Pyjama cover

INTENDED USE: Light fastening to prevent falls and access to the diaper, probes, tubes…

BASIC UDI-DI: 843603535BFXKB

SIZE CHART

Ref.	Fabric	Size	Measures (cm)
BFX8000 BFX8050	Warm Cold	One size	For single beds 200 x 90

MODEL

BFX8000:
BFX8050:

COLOUR

Light blue / white
White

COMPOSITION

40% CO - 60% PES
50% CO - 50% PES
LATEX FREE

LAUNDRY INSTRUCTIONS **BFX8000:**



- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Bodyfix - Pyjama cover

LAUNDRY INSTRUCTIONS **BFX8050:**



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Ironing at low temperature
- Do not dry clean



The product BFX8050 can be machine washed up to 90ºC. However, it will cause a shrinkage of the tissue and may accelerate the degradation of the product.

Shrinkage at 90ºC (First wash):
Weft >15% / Warp >8%

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE