BODYSYSTEM

ENGLISH




Bodysystem pyjamas *have been designed to prevent access to the diaper by the patient. Moreover, due to the easy access without removing the pyjama, sanitization tasks are quick and easy for the caregiver.*

Made of comfortable and breathable fabric, zippered inside the legs (ankle to ankle), on the back or both areas at once, depending on the model.

Available in various colors, white, blue, or pink.

PYJAMA

FEATURES



4 Sizes



2 models:
- Leg zip
- Back/Leg zip



Short / Long sleeves



Soft and breathable Interlock fabric
White colour
Blue colour
Pink colour

SPECIFICATIONS

PRODUCT FAMILY: Bodysystem

PRODUCT NAME: Pyjama

INTENDED USE: Prevents access to the diaper by the patient and facilitate hygiene tasks for caregivers

BASIC UDI-DI: 843603535BSYLL

SIZE CHART

Leg zip							
LONG SLEEVES	SHORT SLEEVES	LONG SLEEVES	SHORT SLEEVES	LONG SLEEVES	SHORT SLEEVES	Size	Waist contour (cm)
BSY8100	BSY8105	BSY8100A	BSY8105A	BSY8100R	BSY8105R	S	100
BSY8110	BSY8115	BSY8110A	BSY8115A	BSY8110R	BSY8115R	M	104
BSY8120	BSY8125	BSY8120A	BSY8125A	BSY8120R	BSY8125R	L	108
BSY8130	BSY8135	BSY8130A	BSY8135A	BSY8130R	BSY8135R	XL	112

Back / Leg zip							
LONG SLEEVES	SHORT SLEEVES	LONG SLEEVES	SHORT SLEEVES	LONG SLEEVES	SHORT SLEEVES	Size	Waist contour (cm)
BSY8300	BSY8305	BSY8300A	BSY8305A	BSY8300R	BSY8305R	S	100
BSY8310	BSY8315	BSY8310A	BSY8315A	BSY8310R	BSY8315R	M	104
BSY8320	BSY8325	BSY8320A	BSY8325A	BSY8320R	BSY8325R	L	108
BSY8330	BSY8335	BSY8330A	BSY8335A	BSY8330R	BSY8335R	XL	112

COLOUR

- White
- Blue
- Pink

COMPOSITION

50% Cotton - 50% Polyester
Polyamide zip
LATEX FREE

PACKAGING



- 1 - Bodysystem - Pyjama

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Ironing at low temperature
- Do not dry clean



Pyjamas can be machine washed up to 90ºC. However, it will cause a shrinkage of the tissue and may accelerate the degradation of the product.

Shrinkage at 90ºC (First wash):
Weft >15% / Warp >8%

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE