
MEDICOVERS

ENGLISH





Medicovers *hygienic cover for mattress protection. Available with and without zipper.*

Made of polyurethane coated polyester fabric. Impermeable, breathable, fireproof, anti-bacterial, anti-fungal and mite barrier construction.
White nylon zip fastener (Only zip products)

MATTRESS COVER

FEATURES

- Oeko-Tex® Standard 100 - Class I
- Fireproof
- Impermeable
- Breathable
- Washable 95ºC
- Autoclave Sterilization

SPECIFICATIONS

PRODUCT FAMILY: MEDICOVERS

PRODUCT NAME: Mattress cover

INTENDED USE: Mattress protection

BASIC UDI-DI: 843603535FPULP

SIZE CHART

Code		Measures (cm)		
Without zip	**With zip**	**Length**	**Width**	**Height**
FPU5000SAN	FPU5200SAN	190	90	15
FPU5010SAN	FPU5210SAN	190	90	18
FPU5020SAN	FPU5220SAN	190	90	20

COLOUR

Light blue

COMPOSITION

38% Polyester - 62% Polyurethane
LATEX FREE

TECHNICAL SPECS

- **Certificate:** Oeko-Tex® Standard 100 Product class I. Tested for harmful substances.
- **Fireproof specs:** Standard BS 7175:1989 Section three. Individual materials.
- **Impermeability:** 180 cm water column
- **Breathability:** 1000 gr/m2 / 24 hours
- **Shrinkage at 95ºC (first wash):** Weft <4% / Warp <2%
- **Autoclave Sterilization:** 120/130 ºC

LAUNDRY INSTRUCTIONS



- Machine wash 95ºC
- Do not bleach
- Tumble dry low
- Iron low (from reverse only)
- Dry Clean (No Trichloroethylene)

- Do not use products with peroxides, bleach, chlorine, dimethyl-formamide, ketone, ammoniac, derivates, phenols, acids y alkalis.

PACKAGING



- 1 Medicovers - Mattress cover

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE