MEDICOVERS

ENGLISH



Medicovers *hygienic cover for pillow protection. Available with and without zipper.*

Made of polyurethane coated polyester fabric. Impermeable, breathable, fireproof, anti-bacterial, anti-fungal and mite barrier construction. White nylon zip fastener (Only zip products).

PILLOW COVER

FEATURES

- Oeko-Tex® Standard 100 - Class I
- Fireproof
- Impermeable
- Breathable
- Washable 95ºC
- Autoclave Sterilization

SPECIFICATIONS

PRODUCT FAMILY: MEDICOVERS

PRODUCT NAME: Pillow cover

INTENDED USE: Pillow protection

BASIC UDI-DI: 843603535FPULP

SIZE CHART

Code		Measures (cm)	
With zip	**Without zip**	**Length**	**Width**
FPU5300SAN	FPU5305SAN	100	43
FPU5310SAN	FPU5315SAN	90	43
FPU5320SAN	FPU5325SAN	80	43
FPU5330SAN	FPU5335SAN	70	43

COLOUR

Light blue

COMPOSITION

38% Polyester - 62% Polyurethane
LATEX FREE

TECHNICAL SPECS

- **Certificate:** Oeko-Tex® Standard 100 Product class I. Tested for harmful substances.
- **Fireproof specs:** Standard BS 7175:1989 Section three. Individual materials.
- **Impermeability:** 180 cm water column
- **Breathability:** 1000 gr/m2 / 24 hours
- **Shrinkage at 95ºC (first wash):** Weft <4% / Warp <2%
- **Autoclave Sterilization:** 120/130 ºC

LAUNDRY INSTRUCTIONS



- Machine wash 95ºC
- Do not bleach
- Tumble dry low
- Iron low (from reverse only)
- Dry Clean (No Trichloroethylene)

- Do not use products with peroxides, bleach, chlorine, dimethyl-formamide, ketone, ammoniac, derivates, phenols, acids y alkalis.

PACKAGING



- 1 Medicovers - Pillow cover

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE