MEDICOVERS

ENGLISH



Medicovers bedrail protector *avoids injuries to the patient by protecting them from blows against the handrail and prevent from inserting the limbs between the lateral bars. With fast assembly and disassembly thanks to the 4 fastening (140 cm / 190 cm) and 2 fastening (100 cm) straps with clip closure .*

Bedrail protector cover made of polyurethane coated polyester fabric. Impermeable, breathable, fireproof, anti-bacterial, anti-fungal and mite barrier construction. Interior made of 25mm width polyurethane foam. Removable for washing. Polypropylene straps and adjustable clip closure. White nylon zipper.



BEDRAIL PROTECTOR

FEATURES

- Oeko-Tex® Standard 100 - Class I
- Fireproof
- Impermeable
- Breathable
- Washable 95ºC
- Autoclave Sterilization

SPECIFICATIONS

PRODUCT FAMILY: MEDICOVERS

PRODUCT NAME: Bedrail protector

INTENDED USE: Prevention of contusions

PRODUCT NAME: Bedrail protector

SIZE CHART

Code	Measures (cm)		
	Length	Width	Thickness
FPU5710 SAN	100	34	3
FPU5740SAN	140	34	3
FPU5750 SAN	190	34	3

COLOUR

- Dark blue

COMPOSITION

Cover: 38% Polyester - 62% Polyurethane
Foam: Polyurethane
LATEX FREE

TECHNICAL SPECS

- **Certificate:** Oeko-Tex® Standard 100 Product class I. Tested for harmful substances.
- **Fireproof specs:** Standard BS 7175:1989 Section three. Individual materials.
- **Impermeability:** 180 cm water column
- **Breathability:** 1000 gr/m2 / 24 hours
- **Shrinkage at 95ºC (first wash):** Weft <4% / Warp <2%
- **Autoclave Sterilization:** 120/130 ºC

LAUNDRY INSTRUCTIONS



- Machine wash 95ºC
- Do not bleach
- Tumble dry low
- Iron low (from reverse only)
- Dry Clean (No Trichloroethylene)

- Do not use products with peroxides, bleach, chlorine, dimethylformamide, ketone, ammoniac, derivates, phenols, acids y alkalis.

PACKAGING



- 1 - Medicovers - Bedrail protector (140/190cm)
- 2 - Medicovers - Bedrail protector (100cm)

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice