SALVACLIP Comfort

ENGLISH



Salvaclip Comfort vest for wheelchair *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. It adapts to the patient's body ergonomically, providing a great feeling of comfort.*

Its ergonomic design attaches to the wheelchair via high-strength straps that fasten with clips. It features reinforced construction with high-tenacity loops at the points with the highest load demand.


Made of tridimensional polyester material, it offers high breathability, support, and a hydrophobic surface that prevents stains from penetrating the fabric and facilitates cleaning. Reinforced polyester straps are 45mm wide.

WHEELCHAIR VEST

FEATURES



Acetal side release plastic buckles




Breathable / Comfortable 3D mesh fabric structure

Available in 3 sizes
S - M - L



Wheelchair identification label



Anti stain liquid repellent



Reinforced polyester straps



Back attachment

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Comfort

PRODUCT NAME: Wheelchair vest

INTENDED USE: To help positioning of seated patients

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures (cm) *
		Contour Min./Max.
SCC4104	S	90 / 115
SCC4105	M	100 / 135
SCC4106	L	115 / 150

* The contour measurements are total (wheelchair + patient)

COLOUR

- Blue
- Black straps

COMPOSITION

Vest: 95% Polyester - 5% Polyurethane
Straps: Polyester
Buckles: Acetal (POM)
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Comfort - Wheelchair vest
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE