SALVACLIP Comfort

ENGLISH



Salvaclip Comfort belt for wheelchair *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability.*

It uses a robust plastic buckles as locking device. This system allows a quick and easy application and can easily adjust the product depending on the size of the patient. It is the best choice for use in patients at risk of slipping and falling off the chair, those who adopt incorrect postures or have uncontrolled movements.


Made using 100% Polyester tridimensional mesh structure, providing high breathability, tear and washing resistance.

WHEELCHAIR BELT

FEATURES



Acetal side release plastic buckles




Breathable / Comfortable 3D mesh fabric structure

Available in 2 sizes
M - L



Wheelchair identification label



Anti stain liquid repellent



Reinforced polyester straps



Back attachment

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Comfort

PRODUCT NAME: Wheelchair belt

INTENDED USE: To help positioning of seated patients

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures (cm) Contour Min./Max.
SCC4305	M	90 / 135
SCC4306	L	115 / 165

* The contour measurements are total (wheelchair + patient)

COLOUR

- Blue
- Black straps

COMPOSITION

Belt: 95% Polyester - 5% Polyurethane
Straps: Polyester
Buckles: Acetal (POM)
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Comfort - Wheelchair belt
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE