SALVACLIP Safe

ENGLISH



Salvaclip Safe perineal vest for wheelchair *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. It adapts to the patient's body ergonomically, providing a great feeling of comfort. Perineal support prevents the patient from moving downward.*

It uses mechanical button with magnetic key opening as secure locking system. All tapes are placed toward the back of the chair, thus avoiding any manipulation of the closure by unauthorized persons. It is the safest option to keep an agitated patient controlled, always under health personnel supervision.


Made of tridimensional polyester material, it offers high breathability, support, and a hydrophobic surface that prevents stains from penetrating the fabric and facilitates cleaning. Reinforced polyester straps are 45mm wide

PERINEAL WHEELCHAIR VEST

FEATURES



Magnetic release safety lock system




Breathable / Comfortable 3D mesh fabric structure

Available in 3 sizes
S - M - L



Perineal band vest



Anti stain liquid repellent



Reinforced polyester straps



Back attachment



Wheelchair identification label

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Safe

PRODUCT NAME: Perineal wheelchair vest

INTENDED USE: To help positioning of seated patients

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures (cm) *
		Contour Min./Max.
SCF4164	S	90 / 115
SCF4165	M	100 / 135
SCF4166	L	115 / 150

* The contour measurements are total (wheelchair + patient)

COLOUR

- Blue
- Black straps

COMPOSITION

Vest: 95% Polyester - 5% Polyurethane
Straps: 100% Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Safe - Perineal wheelchair vest
- 1 - Safety lock
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE