SALVACLIP Safe

ENGLISH



Salvaclip Safe vest for wheelchair *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. It adapts to the patient's body ergonomically, providing a great feeling of comfort.*

It uses mechanical button with magnetic key opening as secure locking system. All tapes are placed toward the back of the chair, thus avoiding any manipulation of the closure by unauthorized persons. It is the safest option to keep an agitated patient controlled, always under health personnel supervision.


Made of tridimensional polyester material, it offers high breathability, support, and a hydrophobic surface that prevents stains from penetrating the fabric and facilitates cleaning. Reinforced polyester straps are 45mm wide.

WHEELCHAIR VEST

FEATURES



Magnetic release safety lock system




Breathable / Comfortable 3D mesh fabric structure

Available in 3 sizes
S - M - L



Wheelchair identification label



Anti stain liquid repellent



Reinforced polyester straps



Back attachment

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Safe

INTENDED USE: To help positioning of seated patients

PRODUCT NAME: Wheelchair vest

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures (cm) * Contour Min./Max.
SCF4154	S	90 / 115 cm
SCF4155	M	100 / 135 cm
SCF4156	L	115 / 150 cm

* The contour measurements are total (wheelchair + patient)

COLOUR

- Blue
- Black straps

COMPOSITION

Vest: 95% Polyester - 5% Polyurethane
Straps: 100% Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Safe - Wheelchair vest
- 1 - Safety lock
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE