
SALVACLIP Safe

ENGLISH



Extension band for enlarging the length of Salvaclip Safe and attach it to bigger armchairs.

Made using 100% Polyester straps and stainless steel eyelets.

EXTENSION BAND

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Safe

PRODUCT NAME: Extension band

INTENDED USE: Extension band for enlarging the lenght of Salvaclip Safe and attach it to bigger armchairs

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Measures
SCF4754	40 cm
SCF4756	60 cm

COLOUR

- Black

COMPOSITION

Strap: 100% Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 5 uni. Salvaclip Safe - Extension band

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE