SALVACLIP

ENGLISH




Salvaclip Standard classic perineal vest *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. Perineal support prevents the patient from moving downward.*

This vest is suitable for wheelchair and armchair. Straps without buckles, should be fastened by knots.


Made of polyester with water repellent coating and comfortable polyurethane foam. Reinforced polyester straps.

CLASSIC PERINEAL VEST

FEATURES



Anti stain liquid repellent



Reinforced polyester straps



Suitable for wheelchair and armchair

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Standard

INTENDED USE: To help positioning of seated patients

PRODUCT NAME: Classic perineal vest

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Straps length	Perineal straps length
SCP4050	One size	4 units - 125 cm	2 units - 175 cm

COLOUR

- Ext.Blue - Int.Black
- Dark blue straps

COMPOSITION

Vest: Ext: Polyester - Int: Polyurethane
Straps: Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Standard - Classic perineal vest
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE