SALVACLIP

ENGLISH



Salvaclip Standard perineal Velcro® belt *to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. Perineal support prevents the patient from moving downward.*

Side release plastic buckles and hook/loop fasteners allows a quick and easy application and can easily adjust the product depending on the size of the patient. It is the best choice for use in patients at risk of slipping and falling off the chair, those who adopt incorrect postures or have uncontrolled movements.


Made of polyester with water repellent coating and comfortable polyurethane foam. Reinforced polyester straps.

PERINEAL VELCRO® BELT


FEATURES



Side release plastic buckles



Anti stain liquid repellent



Reinforced polyester straps



Suitable for wheelchair and armchair



16cm width waist belt with hook/loop fasteners

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Standard

INTENDED USE: To help positioning of seated patients

PRODUCT NAME: Perineal Velcro® belt

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures (cm) Contour Min./Max.
SCP4275	M	100 / 140
SCP4275L	L	120 / 165

COLOUR

- Ext.Blue - Int.Black
- Dark blue straps

COMPOSITION

Ext: Polyester - Int: Polyurethane
Straps: Polyester
Buckles: Nylon
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Standard - Perineal Velcro® belt
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE