SALVACLIP

ENGLISH




Salvaclip Standard waist belt *intended to prevent falling out from bed or slipping down the armchair for elderly, disoriented and non-aggressive patients. It can also be used for aiding the caregiver while making patient transfers between bed and armchair. It includes features from safety restraints systems, such as; Perineal band, sidebands and articulated joint. Acetal triple-lock buckles and Velcro® straps provide a fast and easy adjustment.*


Made of soft polyurethane foam and coated with anti-stain waterproof treatment polyester fabric on both sides. Reinforced with polyester webbing that increases its durability and offering good dimensional stability when laundering.

BED-ARMCHAIR WAIST BELT

FEATURES



Anti stain liquid repellent



Reinforced polyester straps



Articulated joint for a controlated movement



Safety side bands



Safety double lock buckle



Sewn Perineal band

Available in 3 sizes
S - M - L

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Standard

PRODUCT NAME: Bed-Armchair waist belt

INTENDED USE: To help positioning of seated patients

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	Size	Measures Min./Max. (cm)	
		Waist belt	Sewn perineal
SCP4511P	S	45 - 80	60 - 100
SCP4510P	M	70 - 105	60 - 100
SCP4512P	L	95 - 130	60 - 100

COLOUR

- Exterior Blue
- Interior Black
- Dark blue straps

COMPOSITION

Perineal waist belt:
Exterior: Polyester
Interior: Polyurethane

Straps: Polyester

Buckles: Acetal (POM)
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvaclip Standard - Perineal waist belt
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE