SALVACLIP

ENGLISH



Salvaclip Standard belt with pelvic support *designed to prevent seated patients from sliding off the seat. Especially useful for patients with stability problems due to a functional disability.*

Plastic buckles for quick and easy adjustment to the size of the patient.


Made of polyester with water repellent coating and comfortable polyurethane foam. Reinforced polyester straps.

BELT WITH PELVIC SUPPORT


FEATURES



Side release plastic buckles



Anti stain liquid repellent



Reinforced polyester straps



Suitable for wheelchair and armchair

SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Standard

PRODUCT NAME: Belt with pelvic support

INTENDED USE: to help positioning of seated patients

BASIC UDI-DI: 843603535SCPM7

SIZE CHART

Ref.	MEASURES (cm)
SCP4650	PELVIC: Straps length Min. 8 / Max. 125 BELT: Contour Min. 90 / Max. 165

COLOUR

- Ext.Blue - Int.Black
- Dark blue straps

COMPOSITION

Ext: Polyester - Int: Polyurethane
Straps: Polyester
Buckles: Nylon
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Clip belt
- 1 - Pelvic support
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE