SLING ECO

ENGLISH



Sling Eco *is an economical, soft-touch shoulder immobilizer for patients of all ages. It supports and restricts shoulder movement and can be applied to both arms.*

Made of a special padded velour fabric that allows Velcro® to adhere across its entire surface using straps of the same material, which serve to adjust and immobilize the arm.
Sandwich construction with 3 mm foam, non-woven white fabric on the inside, and soft beige velour on the outside. Ultrasonically welded for a seamless construction. White plastic rings 30 mm wide. Beige Velcro® 38 mm wide.

SLING ECO SHOULDER IMMOBILIZER

SPECIFICATIONS

PRODUCT FAMILY: SLING

PRODUCT NAME: SLING ECO shoulder immobilizer

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Measures (cm)	
		Sling	Straps
SOF0009	XS	30 x 12	65 x 4,5 (x2)
SOF0010	S	35 x 12	75 x 4,5 (x2)
SOF0011	M	45 x 14	95 x 4,5 (x2)
SOF0012	L	50 x 14	115 x 4,5 (x2)
SOF0013	XL	50 x 14	135 x 4,5 (x2)

COLOUR

Beige / White

COMPOSITION

35% Polyamide
35% Polyester
30% Viscose
Polyamide D shaped Rings
LATEX FREE

LAUNDRY INSTRUCTIONS

- Do not wash
- Disposable Product

PACKAGING

CARTON BOX 10 UNITS

- 1 - SLING ECO Shoulder immobilizer
- 1 - Instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE