CHARPA Eco

ENGLISH



Charpa Eco *is an inexpensive arm-shoulder immobilizer for patients of all ages. Supports the arm and limits the shoulder movement. Can be placed in both arms.*

Made of non-woven arm support and comfortable padded straps made of (Velour/Foam/Viscose) fabric ultrasonically welded which provides a seamless construction. Its soft touch velour permits the adhesion of hook fasteners over the entire surface, allowing a complete adaptability to patients of diverse complexities.

CHARPA ECO SHOULDER IMMOBILISER

FEATURES



Hook/loop fastening



Polyamide D shaped Rings



Padded thumb support

SPECIFICATIONS

PRODUCT FAMILY: SOFT GOODS

PRODUCT NAME: CHARPA ECO Shoulder immobiliser

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Measures (cm)	
		Body	Straps
SOF0019	XS	33 x 13	68 x 5 (x2)
SOF0020	S	41 x 16	88 x 5 (x2)
SOF0021	M	46 x 19	108 x 5 (x2)
SOF0022	L	50 x 20	118 x 5 (x2)
SOF0023	XL	53 x 22	128 x 5 (x2)

COLOUR

Dark blue / Black

COMPOSITION

30% Polyester
20% Polyamide
50% Polypropylene
Polyamide D shaped Rings
LATEX FREE

LAUNDRY INSTRUCTIONS

- Do not wash
- Disposable product

PACKAGING

CARTON BOX 10 UNITS

- 1 - CHARPA ECO Shoulder immobiliser in plastic bag (individually packed)

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE