SLING roll

ENGLISH



Sling Roll *is a universal system in roll format for securing and immobilizing with Velcro® adjustment. Its construction of soft and padded material provides comfortable and durable support. The roll format offers wide versatility of use, making it suitable as a sling (cast support), arm-shoulder immobilizer, strap for monitoring equipment, securing catheters, and more. Versatile and adaptable for patients of all ages. Easy, quick, and simple to apply. Seam-free surface, made with ultrasonic welding.*

Made of 3 mm thick navy blue velour sandwich fabric, allowing Velcro® to adhere across its entire surface. Navy Velcro® 38 mm wide.

SLING ROLL

SPECIFICATIONS

PRODUCT FAMILY: SLING

PRODUCT NAME: SLING Roll

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Straps	Fasteners
SOF0050	50 meters (3roll x 5cm width)	55 hook fasteners (38mm x 10cm)

COLOUR

- Dark blue

COMPOSITION

35% Polyamide
35% Polyether
30% Polypropylene
LATEX FREE

LAUNDRY INSTRUCTIONS

- Do not wash
- Disposable product

PACKAGING

CARTON BOX 50 METERS

- 50 meter (3 rolls)
- 55 hook fasteners (38mm x 10cm)

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE