SLING

ENGLISH



Sling *- Universal System for Shoulder Immobilization. Its solid construction makes it a safe and comfortable support for the patient. Versatile and adaptable, it can be used by adult patients of all ages, on both the right and left arm. Easy and quick to apply. Seam-free surface with ultrasonic welding.*

Made of a special velour fabric, padded beige color that allows Velcro® to adhere across its entire surface using straps of the same material, which serve to adjust and immobilize the arm. Sandwich construction with 4 mm beige foam, with a weight of 200 gr/m².
White plastic rings 30 mm wide. Beige Velcro® 38 mm wide.

SLING SHOULDER IMMOBILISER

SPECIFICATIONS

PRODUCT FAMILY: SOFT GOODS

PRODUCT NAME: SLING Shoulder immobiliser

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Measures (cm)	
		Sling	Straps
SOF0111	S	35 x 13	(x2) 75 x 6
SOF0101	M	45 x 14	(x2) 95 x 6
SOF0121	L	50 x 14	(x2) 115 x 6
SOF0131	XL	50 x 14	(x2) 135 x 6

COLOUR

- Beige

COMPOSITION

72% Polyester
28% Poliamide
Polyamide D shaped rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Sling - Shoulder immobiliser
- 1 - Instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE