CHARPA

ENGLISH



***Charpa** is a vented light arm and shoulder immobilizer for patients of all ages. It supports the arm and limits shoulder movement. It can be placed on either arm (Ambidextrous). Breathable and washable.*

Made with 100% polyester breathable mesh fabric and comfortable padded straps made of 8mm thick polyester foam, with an inner face of viscose fiber and an outer face of soft polyamide velour ultrasonically welded, providing a seamless construction. Its soft-to-the-touch velour allows the Velcro® to adhere to the entire surface, allowing complete adaptability to patients of various complexities. White polyamide D-rings.

CHARPA SHOULDER IMMOBILISER

SPECIFICATIONS

PRODUCT FAMILY: SOFT GOODS

PRODUCT NAME: CHARPA Shoulder immobiliser

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI : 843603535SOFMP

SIZE CHART

Ref.	Size	Measures (cm)	
		Body	Straps
SOF0201	S	35 x 16	68 x 6
SOF0202	M	45 x 20	94 x 6

COLOUR

White / Light blue

COMPOSITION

Body: 100% Polyester
Straps: 65% Polyester
30% Polyamide
5% Viscose
Polyamide D shaped rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - CHARPA shoulder immobiliser (individually packed)

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE