CHARPA Pro

ENGLISH



***Charpa Pro** is a lightweight arm-shoulder immobilizer for patients of all ages. Supports the arm and limits the shoulder movement. Can be placed in both arms. Breathable and washable.*

Made of polyester fabric arm support and comfortable padded straps made of (Velour / Foam / Non woven) fabric ultrasonically welded which provides a seamless construction. Its soft touch velour permits the adhesion of hook fasteners over the entire surface, allowing a complete adaptability to patients of diverse complexities. Reinforced round rings.

CHARPA PRO SHOULDER IMMOBILISER

FEATURES



Washable
Breathable



Hook/loop fastening



Reinforced
Nickel rings



Padded
thumb support

SPECIFICATIONS

PRODUCT FAMILY: SOFT GOODS

PRODUCT NAME: CHARPA PRO Shoulder immobiliser

INTENDED USE: Shoulder immobilization for treatment of shoulder and arm injuries treatment

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Measures (cm)	
		Body	Straps
SOF0279	XS	33 x 13	72 x 5 (x2)
SOF0280	S	41 x 16	92 x 5 (x2)
SOF0281	M	46 x 19	112 x 5 (x2)
SOF0282	L	50 x 20	122 x 5 (x2)
SOF0283	XL	53 x 22	132 x 5 (x2)

COLOUR

◑ Dark blue / Black

COMPOSITION

Body: 100% Polyester
Straps: 45% Poliuretane
30% Polyamide
25% Viscose
Nickel-free rings
LATEX FREE
Non sterile

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - CHARPA PRO Shoulder immobiliser in plastic bag (individually packed)

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE