ORTOCEL Soft

ENGLISH



OrtocelSoft Knee Immobilizer *recommended for knee injuries, post-operation and post-arthroscopy treatments, immobilization for arthritis, sprains, muscle strain treatment and after premature removal of plaster bandages.*


Made of a three-panelled padded 9mm foam structure that provides an easy adaptation to legs of different diameters, even if the difference between the thigh and the calf are significant. Side panels can be adjusted to the main panel by hook/loop fasteners, easy and quick to apply. Aluminum splints can be placed on each of the panels for temporary immobilization of the knee in extension.

KNEE IMMOBILIZER

FEATURES



Fits both legs



Hook/Loop fastening system



Removable preformed aluminum splints



Adjustable three-panelled structure

SPECIFICATIONS

PRODUCT FAMILY: ORTOCELSOFT

PRODUCT NAME: Knee immobilizer

INTENDED USE: Immobilization for trauma, postoperative and post-arthroscopy treatment, arthritis treatment, dislocations, muscle injuries and after early removal of plaster.

BASIC UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Length (cm)
SOF0400	S	40
SOF0425	M	50
SOF0450	L	60
SOF0475	XL	70

COLOUR

Blue/white

COMPOSITION

55% Polyether 10% Cotton
25% Polyamide 10% Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - OrtocelSoft - Knee immobilizer

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE