ORTOBAND

ENGLISH



Ortoband abdominal binder; *anatomical binder for recovery and prevention of back injuries, post-operative bandaging, usable for pelvic retractions.*

It consists of three compact elastic bands 24 cm wide with a hook/loop closure system.


Available in 5 sizes, differentiated by color code.

ABDOMINAL BINDER

FEATURES



Triple compact elastic band



Hook/Loop fastening system



24cm back height



5 sizes available colour code

SPECIFICATIONS

PRODUCT FAMILY: ORTOBAND

PRODUCT NAME: Abdominal binder

INTENDED USE: Indicated for post-operative recovery and prevention of back injuries

Basic UDI-DI: 843603535SOFMP

SIZE CHART

Ref.	Size	Colour Code	Measures (cm)		
			Length	Stretched contour	
				Min.	Max.
SOF0701	S		86	70	95
SOF0702	M		96	80	105
SOF0703	L		106	90	125
SOF0704	XL		116	100	145
SOF0705	XXL		126	110	165

COLOUR

○ White

COMPOSITION

76%Polyester
17%Elastomer
7%Polyamide

LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 30ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Ortoband - Abdominal binder

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE