SALVACEL

ENGLISH



SVC2400



SVC2420

Limb holder or immobilization system for non-agitated patients by SALVACEL. *The padded and breathable straps provide greater comfort and reduce pressure. Quick and easy installation with Velcro® closures.*

Made of soft-touch padded material (274g/m2). 100% viscose inner coating, breathable 8mm polyester foam, and outer layer of smooth light blue velour fabric, 100% polyamide, adhering to Velcro®. Polyester restraint straps with a breaking strength of up to 160Kg - 1560N, and cotton straps for the SVC2590/SVC2580 model. Quick closure system with Velcro®.



LIMB HOLDER SYSTEM

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL

INTENDED USE: Preventive soft limb holder for non-agitated patients

PRODUCT NAME: Limb holder system

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)	Straps length
SVC2400	Adult	30 x 8	110 x 4
SVC2420	Adult	30 x 8	120 x 4 (2u)
SVC2580	Neonate	15 x 4	80 x 1,5 (2u)
SVC2590	Pediatric	17 x 4	80 x 1,5 (2u)

COLOUR

- Light blue
- Natural white

COMPOSITION


65% PES - 30% PA - 5% Viscose
Polyester straps
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING





- 1 Pair - Salvacel - Limb holder system

ATTENTION

- Single-patient product

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE