SALVACEL

ENGLISH




SALVACEL Abdominal Belt *is designed to securely hold and control the torso while in bed, chair, or wheelchair, suitable for non-agitated patients. It helps prevent slipping and falling from the bed. Its construction with soft and padded materials provides comfortable support without constricting the patient's abdomen. With its crisscross straps at the back, it allows for easy positional changes.*

Manufactured with breathable and hypoallergenic inner and outer coating made of 8mm thick polyester foam, featuring a viscose fiber inner face and a soft polyamide velour outer face (274g/m2). Polyester straps with a breaking strength of up to 160kg - 1560N. Nickel-free metal rings ensure safety and durability.

WAIST BELT

FEATURES



Comfortable padded waist belt



Postural changes allowed



Bed frame fixation by strong polyéster straps & metal rings

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL

INTENDED USE: Preventive soft limb holder for non-agitated patients

PRODUCT NAME: Waist belt

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)	
		Waist Belt	Straps length
SVC2100	S	46 x 12	145
SVC2200	M	67 x 16	170
SVC2300	L	77 x 16	180
SVC2350	XL	115 x 16	205

COLOUR

- Light blue
- Natural white

COMPOSITION


65% PES - 30% PA - 5% Viscose
Polyester straps
Nickel-free metal rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING





- 1 - Salvacel - Waist belt

ATTENTION

- Single-patient product

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE