SALVACEL

ENGLISH



SALVACEL plastic ring limb holder or immobilization system for non-agitated patients*. The padded and breathable limb straps provide greater comfort and reduce pressure. Quick and easy installation with Velcro® closures. The plastic rings provide additional adjustment for the wrists and ankles.*

Made of soft-touch padded material (274gr/m²). 100% viscose inner coating, breathable 8mm polyester foam, and outer layer of smooth light blue velour fabric 100% polyamide, adhering to Velcro®. Polyester restraint straps with a breaking strength of up to 160Kg - 1560N.
Quick closure system with Velcro®.
D-shaped poliamyde ring.



LIMB HOLDER SYSTEM

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL

PRODUCT NAME: Limb holder system

INTENDED USE: Preventive soft limb holder for non-agitated patients

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)	Straps length
SVC2475	Adult	30 x 8	108 x 4

COLOUR

- Light blue
- Natural white

COMPOSITION

65% PES - 30% PA - 5% Viscose
Polyester straps
Polyamide D shaped rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 Pair - Salvacel - Limb holder system

ATTENTION

- Single-patient product

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE