SALVACEL

ENGLISH



SALVACEL mitts *for the restraint of non-agitated patients. Allows securing the patient's hands to the bed to prevent self-injury, aggression, or manipulation of clinical equipment. Quick and easy installation thanks to its hook-and-loop fastening. High-strength polyester mesh fabric for hand restraint, ensuring both durability and breathability.*

Manufactured with a breathable and hypoallergenic inner and outer coating made of 8mm thick polyester foam, with a viscose fiber inner face and a soft polyamide velour outer face (274g/m2). 3mm PEHD splint. Polyester straps with a breaking strength of up to 160Kg - 1560N. Polyamide 'D'-shaped ring.



MITTENS

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL


INTENDED USE: Preventive soft limb holder for non-agitated patients

PRODUCT NAME: Mittens

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)			
		A	B	C	Straps length
SVC2500	M	28	29	16	124
SVC2510	L	31	32	19	124



COLOUR

- Light blue
- Natural white

COMPOSITION

65% PES - 30% PA - 5% Viscose
PEHD inner splint
Polyester straps
Polyamide D-shaped rings
LATEX FREE

PACKAGING





- 1 - Salvacel - Mittens

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

ATTENTION

- Single-patient product

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE