SALVACEL

ENGLISH



SALVACEL protective padded mitten *for the protection of non-agitated patients' hands. It allows for protecting the patient's hands and preventing self-injury, aggression, or manipulation of clinical equipment. Quick and easy installation thanks to Velcro® closure. High-strength and breathable polyester mesh fabric. The fiber interior conforms to the hand's ergonomics, providing greater comfort to the patient.*



Manufactured with a breathable and hypoallergenic inner and outer coating made of 8mm thick polyester foam, with a viscose fiber inner face and a soft polyamide velour outer face (274g/m2). 100% Polyester breathable net fabric. Polyamide 'D'-shaped ring. Hollow siliconized polyester fiber interior. Quick closure system with Velcro®.

PROTECTIVE PADDED MITTEN

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL


PRODUCT NAME: Protective padded mitten

INTENDED USE: Preventive soft limb holder for non-agitated patients

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)		
		A	B	C
SVC2505F	M	28	29	16
SVC2515F	L	31	32	19



COLOUR

- Light blue
- Natural white

COMPOSITION

Mitten: 65% PES - 30% PA - 5% Viscose
Filling: Silicone hollow polyester fibre
Polyamide D-shaped rings
LATEX FREE

PACKAGING



- 1 - Salvacel - Protective padded mitten

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

ATTENTION

- Single-patient product

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE