SALVACEL

ENGLISH



SALVACEL frontal band *incorporates an absorbent zone to avoid the sweat of the forehead during the operations. An elastic band allows easy and quick placement.*

Manufactured from soft-touch padded material. 100% viscose inner lining, breathable 8mm polyester foam padding, and outer velour fabric adherent to Velcro® made of 100% polyamide. Polyester/elastane elastic strap.

FRONTAL BAND

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL

PRODUCT NAME: Frontal band

INTENDED USE: Absorbing sweat from the forehead during operations

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)
SVC2600	One size	20 x 4

COLOUR

- Light blue
- Natural white

COMPOSITION

65% PES - 30% PA - 5% Viscose
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvacel - Frontal band

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE