
SALVAFIX Slim

ENGLISH



SALVAFIX Slim perineal waist belt *for bed with sewn lateral straps to prevent the patient rote over the axis of the belt and sewn perineal band prevent the patient rotate inside the belt and the risk of suffocation by sliding down.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch edging for a comfortable fit. Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet. The mechanism is locked by inserting the pin through its hole and unlocked by placing the magnetic key on its upper face.

PERINEAL WAIST BELT

FEATURES



Adaptable to 80-100cm beds



Soft to the touch edges



Sewn perineal band



Sewn sidebands

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Slim

PRODUCT NAME: Perineal waist belt

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed, stretcher or chair

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size	Measures (cm)	
		Contour Min. / Max.	Width
SVF3050P	M	68 / 108	8

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 1 - Salvafix Slim - Perineal waist belt
- 3 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE