SALVAFIX Transit

ENGLISH



***SALVAFIX Transit wristbands for stretcher**. Its monoblock design permits an easy application on the stretcher.*

Made of 100% polyester straps with inner reinforcement and stainless steel eyelets for a higher strength. 100% polyester soft to the touch piping for a comfortable fit. Welded stainless steel rings. Comfortable padded wristbands with 5mm foam and 100% polyester soft fabric and hook/loop fasteners allows faster and tight application, increasing system security. Safety locks are made of a Stainless Steel pin and a PC/ABS Casing with an internal mechanism that blocks the pin when it is inserted. Magnetic key is needed for system releasing, made of a PA casing and Neodymium magnet.

STRETCHER WRISTBANDS

FEATURES



2 models available:
For 45-60cm and 60-80cm width stretchers



Soft to the touch edges



Comfortable padded wristbands



Colour identifies limb position:
Brown = GROUND (Lower parts)
Light blue = SKY (Upper parts)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Transit

PRODUCT NAME: Stretcher wristbands

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

		Measures (cm)		
Ref.	**Size**	**Stretcher width**	**Contour Min./Max.**	**Width**
SVF3270	M/L	45 - 60	14 / 25	8
SVF3276	M/L	60 - 80	14 / 25	8
SVF3278	XL	60 - 80	24 / 30	8

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 1 - Salvafix Transit - Stretcher Wristbands
- 3 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE