SALVAFIX Transit

ENGLISH



SALVAFIX Transit ankle bands for stretcher. *Safe restraint system with a secure closing by a magnetic button.*
It is suitable for transport agitated and/or disoriented patients on a stretcher.


It prevents patient displacing as well as any possible injuries or aggressions on and with the head.

Designed for an easy and practical installation on the stretcher, eliminating the fixation difficulties of other traditional systems.


It consists of a set of ankle sewn on a band enables the fixation to the strecher.

STRETCHER ANKLE BANDS


FEATURES



2 models available:
For 45-60cm and 60-80cm width stretchers



Soft to the touch edges



Comfortable padded ankle bands



Colour identifies limb position:
Brown = GROUND (Lower parts)
Light blue = SKY (Upper parts)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Transit

PRODUCT NAME: Stretcher ankle bands

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

		Measures (cm)		
Ref.	**Size**	**Stretcher width**	**Contour Min./Max.**	**Width**
SVF3370	M/L	45 - 60	19 / 30	8
SVF3376	M/L	60 - 80	19 / 30	8
SVF3378	XL	60 - 80	29 / 35	8

COLOUR

- Natural white
- Brown

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 1 - Salvafix Transit - Stretcher Ankle bands
- 3 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE