SALVAFIX Transit

ENGLISH



SALVAFIX Thoracic harness. *Safe restraint system with a secure closing by a magnetic button.*


It is suitable for transport agitated and/or disoriented patients on a stretcher.

It prevents patient displacing as well as any possible self- injuries or aggressions with the head.

Designed for an easy and practical installation on the stretcher, eliminating the fixation difficulties of other traditional systems.

STRETCHER HARNESS

FEATURES



2 models available:
For 45-60cm and 60-80cm width stretchers



Soft to the touch edges



8cm width thorax band

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Transit

PRODUCT NAME: Stretcher harness

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the stretcher

BASIC UDI-DI: 843603535SVFNC

SIZE CHART

Code	Size	Measures (cm)		
		Stretcher	Thoracic band	
		Width	Contour Min./Max.	Width
SVF3470	M/L	45 - 60	65 / 110	8
SVF3476	M/L	60 - 80	74 / 114	8
SVF3478	XL	60 - 80	110 / 145	8

COLOUR

- Natural white
- Light blue

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel square rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90ºC
- Do not bleach
- Tumble dry low temperature
- Do not iron

PACKAGING



- 1 - Salvafix Transit - Stretcher Harness
- 3 - Safety locks
- 1 - Magnetic key
- 1 - Safety and instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE