SALVAPAD Mouton

ENGLISH



SALVAPAD Mouton chair cushions *for bedsore prevention. A simple, inexpensive and effective way of preventing sore damage by pressure or abrasion ulcers.*

Cushions stuffed with flakes of hollow silicone fibre. Soft to the touch Mouton fabric on the upper side and polyurethane coated polyester fabric on the bottom, which provides impermeable, breathable, anti-bacterial, anti-fungal and fireproof features. 15mm white straps for chair attachment.

MOUTON CHAIR CUSHIONS

SPECIFICATIONS

PRODUCT FAMILY: SALVAPAD

PRODUCT NAME: Mouton chair cushion

INTENDED USE: Bedsore prevention

Basic UDI-DI: 843603535SVPNY

SIZE CHART

Ref.	Cushion Shape	Measures (cm)
SVP2431	Square	45 x 45 x 10
SVP2432	Square with hole	45 x 45 x 10 (hole ø18)
SVP2433	Horseshoe	45 x 45 x 10
SVP2434	Round with hole	ø45 x 10 (hole ø18)

COLOUR

- Light blue / White
- White straps

COMPOSITION

Upper face: Polyester
Lower face: 62% PU - 38% PES
Fibre: Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 95ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvapad - Mouton chair cushion

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE