SALVAPAD Fibre

ENGLISH





SALVAPAD V-Cushion *for bed sores (decubitus ulcer) prevention and patient positioning. It provides the patient a comfortable and stable position, relieving pressure and improving rest and quality of the patient's sleep. The siliconized hollow fiber interior offers a high-pressure dispersion capacity, making Salvapad Fibre the most convenient and cost-effective choice.*

Made of soft to the touch Cotton/Polyester fabric filled with polyester fiber.

V-CUSHION

SPECIFICATIONS

PRODUCT FAMILY: SALVAPAD Fibre


PRODUCT NAME: V-Cushion

INTENDED USE: Bedsore prevention

BASIC UDI-DI: 843603535SVPNY

SIZE CHART

Code	Size	Measures (cm)		
		A	B	C
SVP2452	One size	70	70	36



COLOUR

- Blue

COMPOSITION

Exterior: 65% Polyester - 35% Cotton
Interior: 100% Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry normal
- Do not iron
- Do not dry clean

ATTENTION

- Vacuum packed: Keep the product out of the bag 24h before its use.

PACKAGING



- 1 - Salvapad Fibre - V-Cushion

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE