SALVAPAD Fibre

ENGLISH





SALVAPAD Fibre chair cushion *for bedsores (decubitus ulcer) prevention. It provides the patient a comfortable and stable position, relieving pressure and improving the rest. Internal cushion with siliconized hollow fibre offers a high-pressure dispersion capacity.*

Removable poliurethane coated cover for an easy maintenance and cleaning. Removable cover made with poliurethane coated polyester fabric on the top and polyester/cotton fabric on the bottom. Hook/loop fastening straps for fixing the cushion to the chair.

CHAIR CUSHION

SPECIFICATIONS

PRODUCT FAMILY: SALVAPAD Fibre

PRODUCT NAME: Chair Cushion

INTENDED USE: Bedsore prevention

BASIC UDI-DI: 843603535SVPNY

SIZE CHART

Code	Size	Measures (cm)
SVP2454	One size	46 x 40 x 8

COLOUR

- Blue

COMPOSITION

Cover:
Top: 51% Polyurethane - 49% Polyester
Bottom: 65% Polyester - 35% Cotton

Cushion:
65% Polyester - 35% Cotton
Interior: 100% Polyester

LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry normal
- Do not iron
- Do not dry clean

ATTENTION

- Vacuum packed: Keep the product out of the bag 24h before its use.

PACKAGING



- 1 - Salvapad Fibre - Chair cushion

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE