SALVAPAD Fibre

ENGLISH



SALVAPAD Fibre wheelchair cushion kit. *Set of 4 cushions for bedsores (decubitus ulcer) prevention. Complete cushion kit designed to fit into wheelchairs, including 1 backrest cushion, 1 seat cushion and 2 armrests cushions. Cushions filled with polyester fiber made of impermeable and breathable fabric on the upper face and comfortable soft to the touch fabric on lower face. It provides the patient, a comfortable and stable position, relieving pressure and improving the rest.*

Hook/loop fastening straps for fixing the cushion to the chair. Zip cover for an easy cleaning and maintenance.

WHEELCHAIR BEDSORE PREVENTION KIT

SPECIFICATIONS

PRODUCT FAMILY: SALVAPAD Fibre

INTENDED USE: Bedsore prevention

PRODUCT NAME: Wheelchair bedsore prevention kit

BASIC UDI-DI: 843603535SVPNY

SIZE CHART

Code	Qty	Cushion	Measures (cm)
SVP2455	1	Square	46 x 40 x 8
	1	Backrest	46 x 40 x 8
	2	Armrests	24 x 29 x 8

COLOUR

- Blue

COMPOSITION

Cover:
Top: 51% Polyurethane - 49% Polyester
Bottom: 65% Polyester - 35% Cotton

Cushion, armrests and backrests:
65% Polyester - 35% Cotton
Interior: 100% Polyester

LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tumble dry normal
- Do not iron
- Do not dry clean

ATTENTION

- Vacuum packed: Keep the product out of the bag 24h before its use.

PACKAGING



- 1 - Salvapad Fibre - Bedsore prevention kit:
  - 1 - Square cushion
  - 1 - Backrest cushion
  - 2 - Armrest cushion

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE