
SALVAPAD Fibre

ENGLISH



SALVAPAD Fibre Cylindrical cushion *provides the patient, a comfortable and stable position, relieving pressure and improving the rest. Internal cushion with siliconized hollow fibre offers a high-pressure dispersion capacity.*

Made of polyurethane coated polyester fabric. Impermeable, breathable, fireproof, anti-bacterial, anti-fungal and mite barrier construction.

CYLINDRICAL CUSHION

FEATURES

- Oeko-Tex® Standard 100 - Class I
- Fireproof
- Impermeable
- Breathable
- Washable 95ºC
- Autoclave Sterilization

SPECIFICATIONS

PRODUCT FAMILY: SALVAPAD Fibre

PRODUCT NAME: Cylindrical cushion

INTENDED USE: Ensures the patient is in a comfortable and stable position

Basic UDI-DI: 843603535SVPNY

SIZE CHART

Code	Measures (cm)	
	Width	Length
SVP2490	Ø20	70

COLOUR

- Light blue

COMPOSITION

Exterior:
38% Polyester - 62% Poliurethane
Interior: 100% Polyester
LATEX FREE

TECHNICAL SPECS

- **Certificate:** Oeko-Tex® Standard 100 Product class I. Tested for harmful substances.
- **Fireproof specs:** Standard BS 7175:1989 Section three. Individual materials.
- **Impermeability:** 180 cm water column
- **Breathability:** 1000 gr/m2 / 24 hours
- **Shrinkage at 95ºC (first wash):** Weft <4% / Warp <2%
- **Autoclave Sterilization:** 120/130 ºC

LAUNDRY INSTRUCTIONS



- Machine wash 95ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

- Do not use products with peroxides, bleach, chlorine, dimethylformamide, ketone, ammoniac, derivates, phenols, acids y alkalis.

ATTENTION

- Product for one patient only

PACKAGING



- 1 Salvapad Fibre - Cylindrical cushion

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice