SALVAFAST

ENGLISH



Salvafast wrist bands *are a fast and easy limb patient restraint. Made in comfortable neoprene and double Hook/Loop fastening system that allows a quick application and release. It prevents the patient from falling out of bed or self-injuring and is appropriate for emergency services, neurology and ICU, where the time factor is crucial. Fixing to the bed by secure closure system to keep the system permanently attached on the bed frame.*

2 options available:

SVT3720 *- Single wristbands: 2 individual pairs.*
SVT3725 *- Monoblock wristbands: Right and left hand sewn in a bed band.*

WRIST BANDS

FEATURES



Safety lock & Magnetic key



Quick attachment Double hook/loop fastening system



Reinforced POM buckles



Permanent bed installation

SPECIFICATIONS

PRODUCT FAMILY: SALVAFAST

PRODUCT NAME: Wristbands

INTENDED USE: Fast restraint of disturbed, violent, aggressive or disoriented patients

Basic UDI-DI: 843603535SVTP8

SIZE CHART

Ref.	Type	Measures (cm)
SVT3720	Single	30 x 6,5
SVT3725	Monoblock	30 x 6,5

COLOUR

Black / Blue

COMPOSITION

80% Polyester - 20% Neoprene
Stainless steel eyelets
Side release POM buckles
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING

- **SVT3720 - Single wristbands**



  - 2 - Salvafast single wristbands
  - 2 - Safety locks
  - 1 - Magnetic key
  - 1 - Instructions manual

- **SVT3720 - Monoblock wristbands**



  - 1 - Salvafast monoblock wristbands
  - 2 - Safety locks
  - 1 - Magnetic key
  - 1 - Instructions manual

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE