SALVAFAST

ENGLISH



Salvafast wrist bands *are a fast and easy limb patient restraint. Made in comfortable neoprene and double Hook/Loop fastening system that allows a quick application and release. It prevents the patient from falling out of bed or self-injuring and is appropriate for emergency services, neurology and ICU, where the time factor is crucial. Fixing to the bed to keep the system permanently attached on the bed frame.*

Single wristbands: 2 individual pairs.

WRIST BANDS

FEATURES



Cushioned soft neoprene



Quick attachment Double hook/loop fastening system



Easy installation with buckle fastener

SPECIFICATIONS

PRODUCT FAMILY: SALVAFAST

PRODUCT NAME: Wristbands

INTENDED USE: Fast restraint of disturbed, violent, aggressive or disoriented patients

Basic UDI-DI: 843603535SVTP8

SIZE CHART

Ref.	Type	Measures (cm)
SVT3722	Single	30 x 6,5

COLOUR

Black / Blue

COMPOSITION

80% Polyester - 20% Neoprene
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 2 - Salvafast single wristbands
- 1 - Instructions manual

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice